

# Instructions of Using DinoXcope Debug Version

This version gives a general description of detection for bug found on normal version.

## Steps

1. To start the debug process, download the DinoXcope 198Debug from <http://www.dino-lite.com/support/DinoXcope198Debug.zip>
2. Connect a Dino-Lite and launch the debug version. Perform where the same issue you have encountered .Then disconnect the Dino-Lite while the application is running. DinoXcope will capture data and display the debug data in the preview comment window.
3. Click in the preview comment window, select all (command + A), copy (command + C) and paste (command + V) in email.

